

产品名称: TR150 Rabbit Polyclonal Antibody
产品货号: APRab19174



产品概述 (Summary)

产品名称 (Production Name)	TR150 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THRAP3 TRAP150
别名 (Alternative Names)	
基因 ID (Gene ID)	9967.0
蛋白 ID (SwissProt ID)	Q9Y2W1.Synthesized peptide derived from human protein . at AA range: 110-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	105kDa

研究背景 (Background)

function:Plays a role in transcriptional coactivation.,PTM:Phosphorylated upon DNA damage, probably by ATM or

产品名称: TR150 Rabbit Polyclonal Antibody
产品货号: APRab19174




ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Subunit of the large multiprotein complex TRAP.,tissue specificity:Ubiquitous.,function:Plays a role in transcriptional coactivation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Subunit of the large multiprotein complex TRAP.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

Signal Transduction; Signaling Pathway; Nuclear Signaling; Nuclear Hormone Receptors; Co-activators/co-repressors; Epigenetics and Nuclear Signaling; Nuclear Signaling Pathways; Nuclear Receptors; Transcription; Mediator Complex

图片 (Image Data)



Western Blot analysis of HeLa lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

注意事项 (Note)

For research use only .