

产品名称: TRIP13 Rabbit Polyclonal Antibody
产品货号: APRab19278



产品概述 (Summary)

产品名称 (Production Name)	TRIP13 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRIP13 TRIP13; PCH2; Pachytene checkpoint protein 2 homolog; Human papillomavirus type 16 E1 protein-binding protein; 16E1-BP; HPV16 E1 protein-binding protein; Thyroid hormone receptor interactor 13; Thyroid receptor-interacting protein 13; TR-in
别名 (Alternative Names)	
基因 ID (Gene ID)	9319.0
蛋白 ID (SwissProt ID)	Q15645. The antiserum was produced against synthesized peptide derived from human TRIP13. AA range: 383-432

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
----------------------	--

产品名称: TRIP13 Rabbit Polyclonal Antibody
产品货号: AP Rab19278



蛋白分子量 (Molecular Weight) 48kDa


研究背景 (Background)

This gene encodes a protein that interacts with thyroid hormone receptors, also known as hormone-dependent transcription factors. The gene product interacts specifically with the ligand binding domain. This gene is one of several that may play a role in early-stage non-small cell lung cancer. [provided by RefSeq, Oct 2009], similarity: Belongs to the AAA ATPase family, subunit: Specifically interacts with the ligand binding domain of the thyroid receptor (TR). This interaction does not require the presence of thyroid hormone for its interaction. Interacts with HPV16 E1.,

研究领域 (Research Area)

Neuroscience; Endocrine system; Thyroid axis

图片 (Image Data)



Western blot analysis of lysates from LOVO and NIH/3T3 cells, using TRIP13 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .