

产品名称: TRMT11 Rabbit Polyclonal Antibody
产品货号: APRab19297



产品概述 (Summary)

产品名称 (Production Name)	TRMT11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRMT11
别名 (Alternative Names)	TRMT11; C6orf75; MDS024; tRNA; guanine(10)-N2)-methyltransferase homolog; tRNA guanosine-2'-O-methyltransferase TRM11 homolog
基因 ID (Gene ID)	60487.0
蛋白 ID (SwissProt ID)	Q7Z4G4.Synthesized peptide derived from the C-terminal region of human TRMT11.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	53kDa

产品名称: TRMT11 Rabbit Polyclonal Antibody
产品货号: APRab19297




研究背景 (Background)

function:Catalytic subunit of an S-adenosyl-L-methionine-dependent tRNA methyltransferase complex that mediates the methylation of the guanosine nucleotide at position 10 (m2G10) in tRNAs.,similarity:Belongs to the methyltransferase superfamily.,function:Catalytic subunit of an S-adenosyl-L-methionine-dependent tRNA methyltransferase complex that mediates the methylation of the guanosine nucleotide at position 10 (m2G10) in tRNAs.,similarity:Belongs to the methyltransferase superfamily.,

研究领域 (Research Area)

Androgen and estrogen metabolism;Histidine metabolism;Tyrosine metabolism;Selenoamino acid metabolism;

图片 (Image Data)



Western blot analysis of HepG2 JK Colo HeLa using TRMT11 antibody.. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .