

产品名称: TSSC1 Rabbit Polyclonal Antibody
产品货号: APRab19373



产品概述 (Summary)

产品名称 (Production Name)	TSSC1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSSC1
别名 (Alternative Names)	TSSC1; Protein TSSC1; Tumor-suppressing STF cDNA 1 protein; Tumor-suppressing subchromosomal transferable fragment candidate gene 1 protein
基因 ID (Gene ID)	7260.0
蛋白 ID (SwissProt ID)	Q53HC9.The antiserum was produced against synthesized peptide derived from human TSSC1. AA range:291-340

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	48kDa

产品名称: TSSC1 Rabbit Polyclonal Antibody
产品货号: AP Rab19373




研究背景 (Background)

This gene has been reported in PMID 9403053 as one of several tumor-suppressing subtransferable fragments located in the imprinted gene domain of 11p15.5, an important tumor-suppressor gene region. Alterations in this region have been associated with the Beckwith-Wiedemann syndrome, Wilms tumor, rhabdomyosarcoma, adrenocortical carcinoma, and lung, ovarian, and breast cancer. Alignment of this gene to genomic sequence data suggests that this gene resides on chromosome 2 rather than chromosome 11. [provided by RefSeq, Dec 2008], similarity: Contains 5 WD repeats.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Cancer susceptibility; Tumor Suppressors

图片 (Image Data)



Western blot analysis of the lysates from HepG2 cells using TSSC1 antibody.

注意事项 (Note)

For research use only .