

产品名称: TSSC4 Rabbit Polyclonal Antibody
产品货号: APRab19375



产品概述 (Summary)

产品名称 (Production Name)	TSSC4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSSC4
别名 (Alternative Names)	
基因 ID (Gene ID)	10078.0
蛋白 ID (SwissProt ID)	Q9Y5U2.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)

This gene is one of several tumor-suppressing subtransferable fragments located in the imprinted gene domain of 11p15.5, an important tumor-suppressor gene region. Alterations in this region have been associated with the Beckwith-Wiedemann


产品名称: TSSC4 Rabbit Polyclonal Antibody
产品货号: APRab19375



syndrome, Wilms tumor, rhabdomyosarcoma, adrenocortical carcinoma, and lung, ovarian, and breast cancer. This gene is located among several imprinted genes; however, this gene, as well as the pan-hematopoietic expression gene (PHEMX), escapes imprinting. This gene may play a role in malignancies and disease that involve this region. Several transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jul 2014],similarity:Belongs to the TSSC4 family.,tissue specificity:Expressed in fetal brain, lung, liver and kidney. Widely expressed in adult tissues.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .