

产品名称: TSSK1 Rabbit Polyclonal Antibody
产品货号: APRab19378



产品概述 (Summary)

产品名称 (Production Name)	TSSK1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TSSK1B SPOGA1 SPOGA4 STK22A
别名 (Alternative Names)	
基因 ID (Gene ID)	83942.0
蛋白 ID (SwissProt ID)	Q9BXA7.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

TSSK1 belongs to a family of serine/threonine kinases highly expressed in testis (Hao et al., 2004 [PubMed 15044604]).
[supplied by OMIM, Mar 2008],catalytic activity:ATP + a protein = ADP + a phosphoprotein.,cofactor:Magnesium.,enzyme

产品名称: TSSK1 Rabbit Polyclonal Antibody
产品货号: APRab19378




regulation: Activated by phosphorylation on Thr-174, potentially by autophosphorylation., function: May be involved in a signaling pathway during male germ cell development or mature sperm

function., PTM: Autophosphorylated., similarity: Belongs to the protein kinase superfamily., similarity: Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family., similarity: Contains 1 protein kinase domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .