

产品名称: UBCE7IP4 Rabbit Polyclonal Antibody
产品货号: APRab19514



产品概述 (Summary)

产品名称 (Production Name)	UBCE7IP4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF144A RNF144A; KIAA0161; RNF144; UBCE7IP4; Probable E3 ubiquitin-protein ligase
别名 (Alternative Names)	RNF144A; RING finger protein 144A; UbcM4-interacting protein 4; Ubiquitin-conjugating enzyme 7-interacting protein 4
基因 ID (Gene ID)	9781.0
蛋白 ID (SwissProt ID)	P50876.The antiserum was produced against synthesized peptide derived from human RNF144A. AA range:151-200

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	30kDa

产品名称: UBCE7IP4 Rabbit Polyclonal Antibody
产品货号: APRab19514




研究背景 (Background)


The protein encoded by this protein contains a RING finger, a motif known to be involved in protein-DNA and protein-protein interactions. The mouse counterpart of this protein has been shown to interact with Ube2l3/UbcM4, which is an ubiquitin-conjugating enzyme involved in embryonic development. [provided by RefSeq, Jul 2008], caution: Lacks the His residue in the RING-type domain 2 that is one of the conserved features of the family., function: E3 ubiquitin-protein ligase which accepts ubiquitin from E2 ubiquitin-conjugating enzymes UBE2L3 and UBE2L6 in the form of a thioester and then directly transfers the ubiquitin to targeted substrates., pathway: Protein modification; protein ubiquitination., similarity: Belongs to the RBR family. RNF144 subfamily., similarity: Contains 1 IBR-type zinc finger., similarity: Contains 2 RING-type zinc fingers., subunit: Interacts with UBE2L3.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from COLO cells, using RNF144A Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using UBCE7IP4 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

产品名称: UBCE7IP4 Rabbit Polyclonal Antibody
产品货号: APRab19514



注意事项 (Note)

For research use only .