

产品名称: Ubiquilin-3 Rabbit Polyclonal Antibody
产品货号: APRab19543



产品概述 (Summary)

产品名称 (Production Name)	Ubiquilin-3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBQLN3
别名 (Alternative Names)	UBQLN3; Ubiquilin-3
基因 ID (Gene ID)	50613.0
蛋白 ID (SwissProt ID)	Q9H347.The antiserum was produced against synthesized peptide derived from human UBQLN3. AA range:271-320

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	70kDa

研究背景 (Background)


产品名称: Ubiquilin-3 Rabbit Polyclonal Antibody
产品货号: APRab19543




Summary: This gene encodes an ubiquitin-like protein (ubiquilin) that shares high degree of similarity with related products in yeast, rat and frog. Ubiquilins contain a N-terminal ubiquitin-like domain and a C-terminal ubiquitin-associated domain. They physically associate with both proteasomes and ubiquitin ligases, and thus thought to functionally link the ubiquitination machinery to the proteasome to affect in vivo protein degradation. This gene is specifically expressed in the testis, and proposed to regulate cell-cycle progression during spermatogenesis. [provided by RefSeq, Jul 2008],similarity:Contains 1 UBA domain.,similarity:Contains 1 ubiquitin-like domain.,tissue specificity:Testis specific.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using UBQLN3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from Jurkat cells using UBQLN3 antibody.

注意事项 (Note)

For research use only .