

产品名称: UBTD2 Rabbit Polyclonal Antibody
产品货号: APRab19588



产品概述 (Summary)

产品名称 (Production Name)	UBTD2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBTD2
别名 (Alternative Names)	UBTD2; DCUBP; SB72; Ubiquitin domain-containing protein 2; Dendritic cell-derived ubiquitin-like protein; DC-UbP; Ubiquitin-like protein SB72
基因 ID (Gene ID)	92181.0
蛋白 ID (SwissProt ID)	Q8WUN7.The antiserum was produced against synthesized peptide derived from human UBTD2. AA range:181-230

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	26kDa

产品名称: UBTD2 Rabbit Polyclonal Antibody
产品货号: APRab19588




研究背景 (Background)


similarity:Contains 1 ubiquitin-like domain.,tissue specificity:Detected in dendritic cells. Highly expressed in tumor cell lines, but not detectable in most tissues.,similarity:Contains 1 ubiquitin-like domain.,tissue specificity:Detected in dendritic cells. Highly expressed in tumor cell lines, but not detectable in most tissues.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of UBTD2 Antibody. The lane on the right is blocked with the UBTD2 peptide.



Western blot analysis of the lysates from HUVEC cells using UBTD2 antibody.

注意事项 (Note)

For research use only .