

产品名称: UEV3 Rabbit Polyclonal Antibody
产品货号: APRab19606



产品概述 (Summary)

产品名称 (Production Name)	UEV3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UEVLD
别名 (Alternative Names)	UEVLD; UEV3; Ubiquitin-conjugating enzyme E2 variant 3; UEV-3; EV and lactate/malate dehydrogenase domain-containing protein
基因 ID (Gene ID)	55293.0
蛋白 ID (SwissProt ID)	Q8IX04.The antiserum was produced against synthesized peptide derived from human UEVLD. AA range:141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	52kDa

产品名称: UEV3 Rabbit Polyclonal Antibody
产品货号: APRab19606




研究背景 (Background)

function: Possible negative regulator of polyubiquitination., similarity: Contains 1 UEV (ubiquitin E2 variant) domain., similarity: In the C-terminal section; belongs to the LDH/MDH superfamily., similarity: In the N-terminal section; belongs to the ubiquitin-conjugating enzyme family. UEV subfamily., subunit: Homodimer., tissue specificity: Colon, colon carcinoma cell lines, normal cervical epithelium, carcinomas of the uterine cervix and peripheral blood leukocytes., function: Possible negative regulator of polyubiquitination., similarity: Contains 1 UEV (ubiquitin E2 variant) domain., similarity: In the C-terminal section; belongs to the LDH/MDH superfamily., similarity: In the N-terminal section; belongs to the ubiquitin-conjugating enzyme family. UEV subfamily., subunit: Homodimer., tissue specificity: Colon, colon carcinoma cell lines, normal cervical epithelium, carcinomas of the uterine cervix and peripheral blood leukocytes.,


研究领域 (Research Area)

Signal Transduction; Protein Phosphorylation; Ser / Thr Kinases; Other Kinases

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using UEVLD Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HepG2 cells, using UEVLD Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: UEV3 Rabbit Polyclonal Antibody
产品货号: APRab19606



注意事项 (Note)

For research use only .