

产品名称: USP38 Rabbit Polyclonal Antibody
产品货号: APRab19678



产品概述 (Summary)

产品名称 (Production Name)	USP38 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP38 USP38; KIAA1891; Ubiquitin carboxyl-terminal hydrolase 38; Deubiquitinating
别名 (Alternative Names)	enzyme 38; HP43.8KD; Ubiquitin thioesterase 38; Ubiquitin-specific-processing protease 38
基因 ID (Gene ID)	84640.0
蛋白 ID (SwissProt ID)	Q8NB14.The antiserum was produced against synthesized peptide derived from human USP38. AA range:551-600

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	120kDa

产品名称: USP38 Rabbit Polyclonal Antibody
产品货号: APRab19678




研究背景 (Background)

catalytic activity: Ubiquitin C-terminal thioester + H₂O = ubiquitin + a thiol., sequence caution: The absence of the residues from Tyr-840 to Ser-871 is not the result of an alternative splicing., similarity: Belongs to the peptidase C19 family., tissue specificity: Highly expressed in skeletal muscle. Expressed in adrenal gland., catalytic activity: Ubiquitin C-terminal thioester + H₂O = ubiquitin + a thiol., sequence caution: The absence of the residues from Tyr-840 to Ser-871 is not the result of an alternative splicing., similarity: Belongs to the peptidase C19 family., tissue specificity: Highly expressed in skeletal muscle. Expressed in adrenal gland.,


研究领域 (Research Area)

Cell Biology; Proteolysis / Ubiquitin; Proteasome / Ubiquitin; Deubiquitination; Epigenetics and Nuclear Signaling; Ubiquitin & Ubiquitin Like Modifiers

图片 (Image Data)




Western blot analysis of USP38 Antibody. The lane on the right is blocked with the USP38 peptide.



Western blot analysis of the lysates from HepG2 cells using USP38 antibody.

产品名称: USP38 Rabbit Polyclonal Antibody
产品货号: APRab19678



Western Blot analysis of various cells using USP38 Polyclonal Antibody diluted at 1 : 1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .