

产品名称: USP42 Rabbit Polyclonal Antibody
产品货号: APRab19680



产品概述 (Summary)

产品名称 (Production Name)	USP42 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	USP42
别名 (Alternative Names)	USP42; Ubiquitin carboxyl-terminal hydrolase 42; Deubiquitinating enzyme 42; Ubiquitin thioesterase 42; Ubiquitin-specific-processing protease 42
基因 ID (Gene ID)	84132.0
蛋白 ID (SwissProt ID)	Q9H9J4.The antiserum was produced against synthesized peptide derived from human USP42. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	150kDa

产品名称: USP42 Rabbit Polyclonal Antibody
产品货号: APRab19680




研究背景 (Background)

catalytic activity: Ubiquitin C-terminal thioester + H₂O = ubiquitin + a thiol, similarity: Belongs to the peptidase C19 family, tissue specificity: Broadly expressed, catalytic activity: Ubiquitin C-terminal thioester + H₂O = ubiquitin + a thiol, similarity: Belongs to the peptidase C19 family, tissue specificity: Broadly expressed,


研究领域 (Research Area)

Cell Biology; Proteolysis / Ubiquitin; Proteasome / Ubiquitin; Deubiquitination; Epigenetics and Nuclear Signaling; Ubiquitin & Ubiquitin Like Modifiers

图片 (Image Data)




Western blot analysis of lysates from HT-29 cells, using USP42 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using USP42 antibody.

产品名称: USP42 Rabbit Polyclonal Antibody
产品货号: APRab19680



Western Blot analysis of various cells using USP42 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .