

产品名称: V1RL2 Rabbit Polyclonal Antibody
产品货号: APRab19698



产品概述 (Summary)

产品名称 (Production Name)	V1RL2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VN1R2
别名 (Alternative Names)	VN1R2; V1RL2; Vomeronasal type-1 receptor 2; G-protein coupled receptor GPCR25; hGPCR25; V1r-like receptor 2
基因 ID (Gene ID)	317701.0
蛋白 ID (SwissProt ID)	Q8NFZ6. The antiserum was produced against synthesized peptide derived from human VN1R2. AA range: 88-137

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	46kDa

产品名称: V1RL2 Rabbit Polyclonal Antibody
产品货号: APRab19698




研究背景 (Background)

function: Putative pheromone receptor., miscellaneous: The chimpanzee, gorilla and orangutan orthologous proteins do not exist, their genes are pseudogenes., online information: No one nose - Issue 61 of August 2005, similarity: Belongs to the G-protein coupled receptor 1 family., function: Putative pheromone receptor., miscellaneous: The chimpanzee, gorilla and orangutan orthologous proteins do not exist, their genes are pseudogenes., online information: No one nose - Issue 61 of August 2005, similarity: Belongs to the G-protein coupled receptor 1 family.,


研究领域 (Research Area)

Neuroscience; Neurotransmission; Receptors / Channels; GPCR; More GPCR; Signal Transduction; Signaling Pathway; G Protein Signaling

图片 (Image Data)



Western blot analysis of lysates from HUVEC, COLO, and MCF-7 cells, using VN1R2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using V1RL2 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .