

产品名称: V1RL5 Rabbit Polyclonal Antibody
产品货号: APRab19700



产品概述 (Summary)

产品名称 (Production Name)	V1RL5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VN1R5
别名 (Alternative Names)	VN1R5; V1RL5; Vomeronasal type-1 receptor 5; G-protein coupled receptor GPCR26; hGPCR26; V1r-like receptor 5
基因 ID (Gene ID)	317705.0
蛋白 ID (SwissProt ID)	Q7Z5H4.The antiserum was produced against synthesized peptide derived from human VN1R5. AA range:44-93

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: V1RL5 Rabbit Polyclonal Antibody
产品货号: APRab19700




研究背景 (Background)

function: Putative pheromone receptor., miscellaneous: The chimpanzee and orangutan orthologous proteins do not exist, their genes are pseudogenes., online information: No one nose - Issue 61 of August 2005, polymorphism: There is an allele VN1R5*2 with a truncating mutation in position 46. It was found with a frequency of about 32%, similarity: Belongs to the G-protein coupled receptor 1 family., function: Putative pheromone receptor., miscellaneous: The chimpanzee and orangutan orthologous proteins do not exist, their genes are pseudogenes., online information: No one nose - Issue 61 of August 2005, polymorphism: There is an allele VN1R5*2 with a truncating mutation in position 46. It was found with a frequency of about 32%, similarity: Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of A549 cells, using VN1R5 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .