

产品名称: Vav2 Rabbit Polyclonal Antibody
产品货号: AP Rab19743



产品概述 (Summary)

产品名称 (Production Name)	Vav2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VAV2
别名 (Alternative Names)	VAV2; Guanine nucleotide exchange factor VAV2; VAV-2
基因 ID (Gene ID)	7410.0
蛋白 ID (SwissProt ID)	P52735.The antiserum was produced against synthesized peptide derived from human VAV2. AA range:108-157

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	101kDa

研究背景 (Background)

产品名称: Vav2 Rabbit Polyclonal Antibody
产品货号: AP Rab19743




VAV2 is the second member of the VAV guanine nucleotide exchange factor family of oncogenes. Unlike VAV1, which is expressed exclusively in hematopoietic cells, VAV2 transcripts were found in most tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2008],function:Guanine nucleotide exchange factor for the Rho family of Ras-related GTPases.,similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,similarity:Contains 1 phorbol-ester/DAG-type zinc finger.,similarity:Contains 1 SH2 domain.,similarity:Contains 2 SH3 domains.,tissue specificity:Widely expressed.,


研究领域 (Research Area)

Regulation of Actin Dynamics; AMPK

图片 (Image Data)




Immunofluorescence analysis of COS7 cells, using VAV2 Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using VAV2 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: Vav2 Rabbit Polyclonal Antibody
产品货号: AP Rab19743



Western blot analysis of lysates from HepG2, HUVEC, and HeLa cells, using VAV2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Vav2 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .