

产品名称: VIP Rabbit Polyclonal Antibody  
产品货号: APRab19803



## 产品概述 (Summary)

|                        |                                |
|------------------------|--------------------------------|
| 产品名称 (Production Name) | VIP Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody     |
| 宿主 (Host)              | Rabbit                         |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA            |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                       |
|------------------------|-------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | VIP                                                                                                   |
| 别名 (Alternative Names) | VIP; VIP peptides                                                                                     |
| 基因 ID (Gene ID)        | 7432.0                                                                                                |
| 蛋白 ID (SwissProt ID)   | P01282.The antiserum was produced against synthesized peptide derived from human VIP. AA range:71-120 |

## 产品应用 (Application)

|                          |                                                                        |
|--------------------------|------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 20kDa                                                                  |

## 研究背景 (Background)

**产品名称: VIP Rabbit Polyclonal Antibody**  
**产品货号: APRab19803**




The protein encoded by this gene belongs to the glucagon family. It stimulates myocardial contractility, causes vasodilation, increases glycogenolysis, lowers arterial blood pressure and relaxes the smooth muscle of trachea, stomach and gall bladder. The protein also acts as an antimicrobial peptide with antibacterial and antifungal activity. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Nov 2014],function:PHM and PHV also cause vasodilation.,function:VIP causes vasodilation, lowers arterial blood pressure, stimulates myocardial contractility, increases glycogenolysis and relaxes the smooth muscle of trachea, stomach and gall bladder.,online information:Vasoactive intestinal peptide entry,similarity:Belongs to the glucagon family.,


### 研究领域 (Research Area)

Neuroscience; Neurotransmitter; Neuropeptides; Hormones

### 图片 (Image Data)



Immunohistochemistry analysis of VIP antibody in paraffin-embedded human testis tissue.



Immunohistochemical analysis of paraffin-embedded Human testis. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody was pre-absorbed by immunogen peptide.

### 注意事项 (Note)

For research use only .

产品名称: VIP Rabbit Polyclonal Antibody  
产品货号: APRab19803

---

