

产品名称: VN1R4 Rabbit Polyclonal Antibody
产品货号: APRab19815



产品概述 (Summary)

产品名称 (Production Name)	VN1R4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VN1R4
别名 (Alternative Names)	VN1R4; V1RL4; Vomeronasal type-1 receptor 4; G-protein coupled receptor GPCR27; hGPCR27; V1r-like receptor 4
基因 ID (Gene ID)	317703.0
蛋白 ID (SwissProt ID)	Q7Z5H5. The antiserum was produced against synthesized peptide derived from human VN1R4. AA range: 88-137

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

产品名称: VN1R4 Rabbit Polyclonal Antibody
产品货号: APRab19815




研究背景 (Background)


function: Putative pheromone receptor., miscellaneous: The chimpanzee and gorilla orthologous proteins do not exist, their genes are pseudogenes., online information: No one nose - Issue 61 of August 2005, polymorphism: Various VN1R4 alleles are known. The sequence shown is that of allele VN1R4*1., similarity: Belongs to the G-protein coupled receptor 1 family., function: Putative pheromone receptor., miscellaneous: The chimpanzee and gorilla orthologous proteins do not exist, their genes are pseudogenes., online information: No one nose - Issue 61 of August 2005, polymorphism: Various VN1R4 alleles are known. The sequence shown is that of allele VN1R4*1., similarity: Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of MCF7 cells, using VN1R4 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from COS7 cells, using VN1R4 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: VN1R4 Rabbit Polyclonal Antibody
产品货号: APRab19815



Western blot analysis of the lysates from COLO205 cells using VN1R4 antibody.

注意事项 (Note)

For research use only .