

产品名称: ZAR1 Rabbit Polyclonal Antibody
产品货号: APRab20029



产品概述 (Summary)

产品名称 (Production Name)	ZAR1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZAR1
别名 (Alternative Names)	ZAR1; Zygote arrest protein 1; Oocyte-specific maternal effect factor
基因 ID (Gene ID)	326340.0
蛋白 ID (SwissProt ID)	Q86SH2.Synthesized peptide derived from the Internal region of human ZAR1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)

This maternal effect gene is oocyte-specific and encodes a protein that is thought to function in the initiation of

产品名称: ZAR1 Rabbit Polyclonal Antibody
产品货号: AP Rab20029




embryogenesis. A similar protein in mouse is required for female fertility. [provided by RefSeq, Jul 2013],function:Essential for female fertility. May play a role in the oocyte-to-embryo transition.,similarity:Belongs to the ZAR1 family.,tissue specificity:Ovary and testis.,

研究领域 (Research Area)

Tags & Cell Markers; Cell Type Markers

图片 (Image Data)



Western Blot analysis of various cells using ZAR1 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .