

产品名称: ZBTB40 Rabbit Polyclonal Antibody
产品货号: APRab20046



产品概述 (Summary)

产品名称 (Production Name)	ZBTB40 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZBTB40
别名 (Alternative Names)	ZBTB40; KIAA0478; Zinc finger and BTB domain-containing protein 40
基因 ID (Gene ID)	9923.0
蛋白 ID (SwissProt ID)	Q9NUA8.The antiserum was produced against synthesized peptide derived from human ZBTB40. AA range:1121-1170

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	140kDa

研究背景 (Background)

产品名称: ZBTB40 Rabbit Polyclonal Antibody
产品货号: APRab20046




function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 BTB (POZ) domain., similarity: Contains 12 C2H2-type zinc fingers., function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 BTB (POZ) domain., similarity: Contains 12 C2H2-type zinc fingers.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Domain Families; Zinc Finger

图片 (Image Data)




Western blot analysis of ZBTB40 Antibody. The lane on the right is blocked with the ZBTB40 peptide.




Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using ZBTB40 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: ZBTB40 Rabbit Polyclonal Antibody
产品货号: APRab20046



Western blot analysis of the lysates from HUVEC cells using ZBTB40 antibody.



Western Blot analysis of various cells using ZBTB40 Polyclonal Antibody. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .