

产品名称: ZC3H1 Rabbit Polyclonal Antibody  
产品货号: APRab20050



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | ZC3H1 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,ELISA                         |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                  |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG                                                                                      |
| 克隆 (Clonality)      | Polyclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化方式 (Purification) | Affinity purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                      |
|------------------------|----------------------------------------------------------------------|
| 基因名 (Gene Name)        | ZFC3H1 CCDC131 KIAA0546 PSRC2                                        |
| 别名 (Alternative Names) |                                                                      |
| 基因 ID (Gene ID)        | 196441.0                                                             |
| 蛋白 ID (SwissProt ID)   | O60293.Synthesized peptide derived from part region of human protein |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 218kDa                               |

## 研究背景 (Background)

PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 5 TPR repeats.,similarity:Contains 6 HAT repeats.,PTM:Phosphorylated upon DNA damage,


产品名称: ZC3H1 Rabbit Polyclonal Antibody  
产品货号: APRab20050



probably by ATM or ATR.,similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 5 TPR repeats.,similarity:Contains 6 HAT repeats.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

## 注意事项 (Note)

For research use only .