

产品名称: ZFP106 Rabbit Polyclonal Antibody
产品货号: APRab20087



产品概述 (Summary)

产品名称 (Production Name)	ZFP106 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZFP106
别名 (Alternative Names)	ZFP106; SH3BP3; ZNF474; Zinc finger protein 106 homolog; Zfp-106; Zinc finger protein 474
基因 ID (Gene ID)	64397.0
蛋白 ID (SwissProt ID)	Q9H2Y7.The antiserum was produced against synthesized peptide derived from human ZFP106. AA range:1833-1882

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	209kDa

产品名称: ZFP106 Rabbit Polyclonal Antibody
产品货号: APRab20087



研究背景 (Background)

PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 2 C2H2-type zinc fingers.,similarity:Contains 6 WD repeats.,subunit:Interacts with TSG118.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Contaminating sequence. Potential poly-A sequence.,similarity:Contains 2 C2H2-type zinc fingers.,similarity:Contains 6 WD repeats.,subunit:Interacts with TSG118.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysate from COS7 cells, using ZFP106 antibody.



Western Blot analysis of various cells using ZFP106 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

产品名称: ZFP106 Rabbit Polyclonal Antibody
产品货号: APRab20087



Western Blot analysis of COS7 cells using ZFP106 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .