

产品名称: ZIP9 Rabbit Polyclonal Antibody
产品货号: AP Rab20114



产品概述 (Summary)

产品名称 (Production Name)	ZIP9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC39A9
别名 (Alternative Names)	SLC39A9; ZIP9; Zinc transporter ZIP9; Solute carrier family 39 member 9; Zrt- and Irt-like protein 9; ZIP-9
基因 ID (Gene ID)	55334.0
蛋白 ID (SwissProt ID)	Q9NUM3.The antiserum was produced against synthesized peptide derived from human SLC39A9. AA range:101-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	33kDa

产品名称: ZIP9 Rabbit Polyclonal Antibody
产品货号: APRab20114




研究背景 (Background)

function: May act as a zinc-influx transporter., similarity: Belongs to the ZIP transporter (TC 2.A.5) family., function: May act as a zinc-influx transporter., similarity: Belongs to the ZIP transporter (TC 2.A.5) family.,


研究领域 (Research Area)

Signal Transduction; Metabolism; Vitamins / Minerals; Pathways and Processes; Cofactors, Vitamins / minerals; Vitamins / minerals


图片 (Image Data)



Western blot analysis of SLC39A9 Antibody. The lane on the right is blocked with the SLC39A9 peptide.



Western blot analysis of the lysates from HeLa cells using SLC39A9 antibody.



Western Blot analysis of various cells using ZIP9 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

产品名称: ZIP9 Rabbit Polyclonal Antibody
产品货号: APRab20114



注意事项 (Note)

For research use only .