

产品名称: ZN-16 Rabbit Polyclonal Antibody
产品货号: APRab20130



产品概述 (Summary)

产品名称 (Production Name)	ZN-16 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF292
别名 (Alternative Names)	ZNF292; KIAA0530; Zinc finger protein 292
基因 ID (Gene ID)	23036.0
蛋白 ID (SwissProt ID)	O60281.The antiserum was produced against synthesized peptide derived from human ZNF292. AA range:1311-1360

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	300kDa

研究背景 (Background)


产品名称: ZN-16 Rabbit Polyclonal Antibody
产品货号: APRab20130



function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 9 C2H2-type zinc fingers.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ZNF292 Antibody. The lane on the right is blocked with the ZNF292 peptide.

注意事项 (Note)

For research use only .