

产品名称: ZNF227 Rabbit Polyclonal Antibody
产品货号: APRab20252



产品概述 (Summary)

产品名称 (Production Name)	ZNF227 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF227
别名 (Alternative Names)	ZNF227; Zinc finger protein 227
基因 ID (Gene ID)	7770.0
蛋白 ID (SwissProt ID)	Q86WZ6.Synthesized peptide derived from the Internal region of human ZNF227.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	92kDa

研究背景 (Background)


产品名称: ZNF227 Rabbit Polyclonal Antibody
产品货号: AP Rab20252




function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 KRAB domain., similarity: Contains 19 C2H2-type zinc fingers., function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 KRAB domain., similarity: Contains 19 C2H2-type zinc fingers.,

研究领域 (Research Area)

图片 (Image Data)




Western Blot analysis of extracts from Jurkat cells, using ZNF227 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .




Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:100

产品名称: ZNF227 Rabbit Polyclonal Antibody
产品货号: APRab20252



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-stomach, antibody was diluted at 1:100

注意事项 (Note)

For research use only .