

产品名称: ZNF232 Rabbit Polyclonal Antibody
产品货号: APRab20254



产品概述 (Summary)

产品名称 (Production Name)	ZNF232 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF232
别名 (Alternative Names)	ZNF232; ZSCAN11; Zinc finger protein 232; Zinc finger and SCAN domain-containing protein 11
基因 ID (Gene ID)	7775.0
蛋白 ID (SwissProt ID)	Q9UNY5. The antiserum was produced against synthesized peptide derived from human ZNF232. AA range: 171-220

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: ZNF232 Rabbit Polyclonal Antibody
产品货号: APRab20254




研究背景 (Background)

function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 SCAN box domain., similarity: Contains 5 C2H2-type zinc fingers., tissue specificity: Ubiquitous. Higher expression seen in the liver, testis and ovary., function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 SCAN box domain., similarity: Contains 5 C2H2-type zinc fingers., tissue specificity: Ubiquitous. Higher expression seen in the liver, testis and ovary.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using ZNF232 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .