

产品名称: ZNF262 Rabbit Polyclonal Antibody
产品货号: APRab20256



产品概述 (Summary)

产品名称 (Production Name)	ZNF262 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZMYM4
别名 (Alternative Names)	ZMYM4; KIAA0425; ZNF262; Zinc finger MYM-type protein 4; Zinc finger protein 262
基因 ID (Gene ID)	9202.0
蛋白 ID (SwissProt ID)	Q5VZL5.The antiserum was produced against synthesized peptide derived from human ZMYM4. AA range:801-850

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	175kDa

产品名称: ZNF262 Rabbit Polyclonal Antibody
产品货号: APRab20256




研究背景 (Background)


miscellaneous:The 3'-UTR region of the mRNA encoding this protein contains a motif called CDIR (for cell death inhibiting RNA) that binds HNRPD/AUF1 and HSPB1/HSP27. It is able to inhibit interferon-gamma induced apoptosis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 9 MYM-type zinc fingers.,tissue specificity:Expressed at higher level in heart, skeletal muscle, kidney and liver.,miscellaneous:The 3'-UTR region of the mRNA encoding this protein contains a motif called CDIR (for cell death inhibiting RNA) that binds HNRPD/AUF1 and HSPB1/HSP27. It is able to inhibit interferon-gamma induced apoptosis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 9 MYM-type zinc fingers.,tissue specificity:Expressed at higher level in heart, skeletal muscle, kidney and liver.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ZMYM4 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

产品名称: ZNF262 Rabbit Polyclonal Antibody
产品货号: APRab20256



注意事项 (Note)

For research use only .