

产品名称: ZNF287 Rabbit Polyclonal Antibody
产品货号: APRab20257



产品概述 (Summary)

产品名称 (Production Name)	ZNF287 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF287
别名 (Alternative Names)	ZNF287; ZKSCAN13; Zinc finger protein 287; Zinc finger protein with KRAB and SCAN domains 13
基因 ID (Gene ID)	57336.0
蛋白 ID (SwissProt ID)	Q9HBT7.The antiserum was produced against synthesized peptide derived from human ZNF287. AA range:241-290

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	90kDa

产品名称: ZNF287 Rabbit Polyclonal Antibody
产品货号: APRab20257




研究背景 (Background)


zinc finger protein 287(ZNF287) Homo sapiens This gene encodes a member of the krueppel family of zinc finger proteins, suggesting a role as a transcription factor. Its specific function has not been determined. This gene is located near the Smith-Magenis syndrome region on chromosome 17. [provided by RefSeq, Jul 2008],function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 14 C2H2-type zinc fingers.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ZNF287 Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from NIH/3T3 cells, using ZNF287 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: ZNF287 Rabbit Polyclonal Antibody
产品货号: APRab20257



Western blot analysis of the lysates from Jurkat cells using ZNF287 antibody.



Western Blot analysis of 3T3 cells using ZNF287 Polyclonal Antibody. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .

注意事项 (Note)

For research use only .