

产品名称: ZNF695 Rabbit Polyclonal Antibody
产品货号: APRab20281



产品概述 (Summary)

产品名称 (Production Name)	ZNF695 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF695
别名 (Alternative Names)	ZNF695; Zinc finger protein 695; Zinc finger protein SBZF3
基因 ID (Gene ID)	57116.0
蛋白 ID (SwissProt ID)	Q8IW36.Synthesized peptide derived from the Internal region of human ZNF695.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	65-70kDa or 20kDa for isoform Q8IW36-1

研究背景 (Background)

产品名称: ZNF695 Rabbit Polyclonal Antibody
产品货号: AP Rab20281




function: May be involved in transcriptional regulation., similarity: Contains 1 KRAB domain., function: May be involved in transcriptional regulation., similarity: Contains 1 KRAB domain.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Domain Families; Zinc Finger

图片 (Image Data)



Western Blot analysis of 293T cells using ZNF695 Polyclonal Antibody diluted at 1: 2000. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .