

产品名称: ZNF771 Rabbit Polyclonal Antibody
产品货号: APRab20286



产品概述 (Summary)

产品名称 (Production Name)	ZNF771 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF771
别名 (Alternative Names)	ZNF771; Zinc finger protein 771; Mesenchymal stem cell protein DSC43
基因 ID (Gene ID)	51333.0
蛋白 ID (SwissProt ID)	Q7L3S4.The antiserum was produced against synthesized peptide derived from human ZNF771. AA range:161-210

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

研究背景 (Background)

产品名称: ZNF771 Rabbit Polyclonal Antibody
产品货号: APRab20286




function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 8 C2H2-type zinc fingers., function: May be involved in transcriptional regulation., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 8 C2H2-type zinc fingers.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Domain Families; Zinc Finger; Transcription Factors


图片 (Image Data)



Western blot analysis of lysates from HepG2 cells, using ZNF771 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using ZNF771 antibody.



产品名称: ZNF771 Rabbit Polyclonal Antibody
产品货号: APRab20286



Western Blot analysis of various cells using ZNF771 Polyclonal Antibody. Secondary antibody was diluted at 1:20000 cells
nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .