

产品名称: ZnT-9 Rabbit Polyclonal Antibody
产品货号: AP Rab20302



产品概述 (Summary)

产品名称 (Production Name)	ZnT-9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC30A9
别名 (Alternative Names)	SLC30A9; C4orf1; HUEL; Zinc transporter 9; ZnT-9; Human embryonic lung protein; HuEL; Solute carrier family 30 member 9
基因 ID (Gene ID)	10463.0
蛋白 ID (SwissProt ID)	Q6PML9.The antiserum was produced against synthesized peptide derived from human SLC30A9. AA range:181-230

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	63kDa

产品名称: ZnT-9 Rabbit Polyclonal Antibody
产品货号: APRab20302




研究背景 (Background)


function: Plays a role in the p160 coactivator signaling pathway that mediates transcriptional activation by nuclear receptors (By similarity). Plays a role in transcriptional activation of Wnt-responsive genes., similarity: Belongs to the cation diffusion facilitator (CDF) transporter (TC 2.A.4) family. SLC30A subfamily., subunit: Interacts with GRIP1, ESR1 and AR., tissue specificity: Ubiquitously expressed in fetal and adult tissues and cancer cell lines., function: Plays a role in the p160 coactivator signaling pathway that mediates transcriptional activation by nuclear receptors (By similarity). Plays a role in transcriptional activation of Wnt-responsive genes., similarity: Belongs to the cation diffusion facilitator (CDF) transporter (TC 2.A.4) family. SLC30A subfamily., subunit: Interacts with GRIP1, ESR1 and AR., tissue specificity: Ubiquitously expressed in fetal and adult tissues and cancer cell lines.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of SLC30A9 Antibody. The lane on the right is blocked with the SLC30A9 peptide.



Western blot analysis of the lysates from K562 cells using SLC30A9 antibody.

产品名称: ZnT-9 Rabbit Polyclonal Antibody
产品货号: APRab20302



Western Blot analysis of HepG2 cells using ZnT-9 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .