

产品名称: ZSCAN22 Rabbit Polyclonal Antibody
产品货号: APRab20322



产品概述 (Summary)

产品名称 (Production Name)	ZSCAN22 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZSCAN22
别名 (Alternative Names)	ZSCAN22; HKR2; ZNF50; Zinc finger and SCAN domain-containing protein 22; Krueppel-related zinc finger protein 2; Protein HKR2; Zinc finger protein 50
基因 ID (Gene ID)	342945.0
蛋白 ID (SwissProt ID)	P10073.The antiserum was produced against synthesized peptide derived from human ZSCAN22. AA range:321-370

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

产品名称: ZSCAN22 Rabbit Polyclonal Antibody
产品货号: APRab20322




研究背景 (Background)


function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 SCAN box domain., similarity: Contains 8 C2H2-type zinc fingers., function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the krueppel C2H2-type zinc-finger protein family., similarity: Contains 1 SCAN box domain., similarity: Contains 8 C2H2-type zinc fingers.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of ZSCAN22 Antibody. The lane on the right is blocked with the ZSCAN22 peptide.



Immunohistochemistry analysis of paraffin-embedded human brain, using ZSCAN22 Antibody. The lane on the right is blocked with the ZSCAN22 peptide.

产品名称: ZSCAN22 Rabbit Polyclonal Antibody
产品货号: APRab20322



Western blot analysis of the lysates from HepG2 cells using ZSCAN22 antibody.

注意事项 (Note)

For research use only .