

产品名称: ZBP1 Rabbit Polyclonal Antibody
产品货号: APRab20347



产品概述 (Summary)

产品名称 (Production Name)	ZBP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZBP1 C20orf183 DLM1
别名 (Alternative Names)	C20orf183 ,DLM1 ,Z-DNA-binding protein 1,DNA-dependent activator of IFN-regulatory factors,Tumor stroma and activated macrophage protein DLM-1
基因 ID (Gene ID)	81030.0
蛋白 ID (SwissProt ID)	Q9H171.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	47kDa

研究背景 (Background)

This gene encodes a Z-DNA binding protein. The encoded protein plays a role in the innate immune response by binding

产品名称: ZBP1 Rabbit Polyclonal Antibody
产品货号: APRab20347




to foreign DNA and inducing type-I interferon production. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011],

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .