

产品名称: CASK Rabbit Polyclonal Antibody
产品货号: APRab20350



产品概述 (Summary)

产品名称 (Production Name)	CASK Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CASK
别名 (Alternative Names)	calcium/calmodulin dependent serine protein kinase; CMG; FGS4; LIN2;TNRC8; CAGH39; CAMGUK; MICPCH; MRXSNA
基因 ID (Gene ID)	8573.0
蛋白 ID (SwissProt ID)	O14936..

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 105 kDa; Observed MW: 105 kDa

研究背景 (Background)

Multidomain scaffolding Mg²⁺-independent protein kinase that catalyzes the phosphotransfer from ATP to proteins such


产品名称: CASK Rabbit Polyclonal Antibody
产品货号: APRab20350




as NRXN1, and plays a role in synaptic transmembrane protein anchoring and ion channel trafficking

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CASK in Human cerebrum lysates using CASK antibody.



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer using CASK antibody.

注意事项 (Note)

For research use only .