

产品名称: SRP68 Rabbit Polyclonal Antibody
产品货号: APRab20351



产品概述 (Summary)

产品名称 (Production Name)	SRP68 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRP68
别名 (Alternative Names)	Signal recognition particle 68 kDa protein
基因 ID (Gene ID)	6730.0
蛋白 ID (SwissProt ID)	Q9UHB9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 71 kDa; Observed MW: 71 kDa

研究背景 (Background)

Component of the signal recognition particle (SRP) complex, a ribonucleoprotein complex that mediates the cotranslational targeting of secretory and membrane proteins to the endoplasmic reticulum (ER).


产品名称: SRP68 Rabbit Polyclonal Antibody
产品货号: APRab20351



研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of SRP68 in 293T, HepG2, K562 lysates using SRP68 antibody

注意事项 (Note)

For research use only .