

产品名称: UDB17 Rabbit Polyclonal Antibody
产品货号: APRab20359



产品概述 (Summary)

产品名称 (Production Name)	UDB17 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	UDB17
别名 (Alternative Names)	UDP-glucuronosyltransferase 2B17, UDPGT 2B17, C19-steroid-specific UDP-glucuronosyltransferase, C19-steroid-specific UDPGT, UGT2B17
基因 ID (Gene ID)	NP_001068.1
蛋白 ID (SwissProt ID)	O75795. Generated from rabbits immunized with synthetic peptide between 119-147 amino acids..

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:500
蛋白分子量 (Molecular Weight)	61kDa

产品名称: UDB17 Rabbit Polyclonal Antibody
产品货号: APRab20359




研究背景 (Background)


UGT2B17 belongs to the family of UDP-glucuronosyltransferases (UGTs; EC 2.4.1.17), enzymes that catalyze the transfer of glucuronic acid from uridine diphosphoglucuronic acid to a variety of substrates, including steroid hormones.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of UDB17 antibody (Center) in Y79 cell line lysates (35ug/lane).



UDB17 Antibody (Center) IHC analysis in formalin fixed and paraffin embedded human Prostate carcinoma followed by peroxidase conjugation of the secondary antibody and DAB staining.

注意事项 (Note)

For research use only .