

产品名称: MOX1 (Meox1) Rabbit Polyclonal Antibody
产品货号: APRab20363



产品概述 (Summary)

产品名称 (Production Name)	MOX1 (Meox1) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS with 0.09% (W/V) sodium azide.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MOX1
别名 (Alternative Names)	Homeobox protein MOX-1, Mesenchyme homeobox 1, MEOX1, MOX1
基因 ID (Gene ID)	NP_001035091.1; NP_004518.1; NP_054705.1
蛋白 ID (SwissProt ID)	P50221..

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	28kDa

研究背景 (Background)


Meox1 is a member of a subfamily of non-clustered, diverged, antennapedia-like homeobox-containing genes. This protein may play a role in the molecular signaling network regulating somite development.

产品名称: MOX1 (Meox1) Rabbit Polyclonal Antibody
产品货号: AP Rab20363



研究领域 (Research Area)

图片 (Image Data)



Western blot analysis in 293 cell line lysates. This demonstrates the Meox1 antibody detected the Meox1 protein.

注意事项 (Note)

For research use only .