

产品名称: IL-17A Rabbit Polyclonal Antibody
产品货号: APRab20367



产品概述 (Summary)

产品名称 (Production Name)	IL-17A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL17A
别名 (Alternative Names)	IL-17A; IL17; CTLA8; IL-17; CTLA-8
基因 ID (Gene ID)	3605.0
蛋白 ID (SwissProt ID)	Q16552.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

The protein encoded by this gene is a proinflammatory cytokine produced by activated T cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and

产品名称: IL-17A Rabbit Polyclonal Antibody
产品货号: APRab20367




cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO). High levels of this cytokine are associated with several chronic inflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human colorectal cancer using IL-17A antibody.

注意事项 (Note)

For research use only .