

产品名称: PBXIP1 Rabbit Polyclonal Antibody
产品货号: APRab20401



产品概述 (Summary)

产品名称 (Production Name)	PBXIP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	pH 7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PBXIP1
别名 (Alternative Names)	antibody; OBPIIa antibody
基因 ID (Gene ID)	57326.0
蛋白 ID (SwissProt ID)	Q96AQ6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 81 kDa; Observed MW: 81 kDa

研究背景 (Background)


Regulator of pre-B-cell leukemia transcription factors (BPXs) function.

产品名称: PBXIP1 Rabbit Polyclonal Antibody
产品货号: APRab20401




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PBXIP1 in PC-3, A549 and Hela lysates using PBXIP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer using PBXIP1 antibody. High-pressure and temperature

Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .