

产品名称: SLUR1 Rabbit Polyclonal Antibody
产品货号: APRab20411



产品概述 (Summary)

产品名称 (Production Name)	SLUR1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
纯化方式 (Purification)	Affinity Purified

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLURP1 ARS
别名 (Alternative Names)	-
基因 ID (Gene ID)	57152.0
蛋白 ID (SwissProt ID)	P55000.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Observed MW:11 kDa

研究背景 (Background)

Caution:It is not certain that ARS and ANUP are identical proteins.,Disease:Defects inSLURP1 are a cause of Mal de Meleda (MDM)[MIM:248300]; also known as keratosispalmoplantaris transgradiens of Siemens. MDM is a rare autosomal recessive

产品名称: SLUR1 Rabbit Polyclonal Antibody
产品货号: APRab20411




skin

disorder, characterized by diffuse transgressive palmoplantar keratoderma with keratotic lesions extending onto the dorsa of the hands and the feet (transgrediens). Patients may have hyperhidrosis. Other features include perioral erythema, lichenoid plaques on the knees and the elbows, and nail abnormalities. Function: Has an antitumor activity. Was found to be a marker of late differentiation of the skin. Implicated in maintaining the physiological and structural integrity of the keratinocyte layers of the skin. Induction: Regulated by retinoic acid, epidermal growth factor and interferon-gamma. Similarity: Contains 1 UPAR/Ly6 domain. Subunit: Homodimer. Tissue specificity: Granulocytes. Expressed in skin. Predominantly expressed in the granular layer of skin, notably the acrosyringium. Identified in several biological fluids such as sweat, saliva, tears, plasma and urine.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4°C overnight

注意事项 (Note)

For research use only.