

产品名称: Wnt3a Rabbit Polyclonal Antibody
产品货号: APRab20422



产品概述 (Summary)

产品名称 (Production Name)	Wnt3a Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS+50%Glycerol+0.05% Proclin 300+0.05%protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	wnt3a MGC119418;MGC119419;MGC119420;ProteinWnt 3a Precursor ;Protein Wnt-3a ;Wingless type MMTV integration site family member 3A ;Wnt3a ;wnt3a ;WNT3A protein ;WNT3A_HUMAN.
别名 (Alternative Names)	
基因 ID (Gene ID)	89780.0
蛋白 ID (SwissProt ID)	P56704.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 39 kDa

研究背景 (Background)

产品名称: Wnt3a Rabbit Polyclonal Antibody
产品货号: APRab20422



The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 96% amino acid identity to mouse Wnt3A protein, and 84% to human WNT3 protein, another WNT gene product. This gene is clustered with

WNT14 gene, another family member, in chromosome 1q42 region. [provided by RefSeq, Jul 2008],


研究领域 (Research Area)

mTOR signaling pathway ;

Wnt signaling pathway ;

Hippo signaling pathway ;

图片 (Image Data)



Various whole cell lysates were separated by 4-20% SDS-PAGE, and the membrane was blotted with anti-Wnt3a antibody.

注意事项 (Note)

For research use only .