

产品名称: PROKR2 Rabbit Polyclonal Antibody
产品货号: APRab20429



产品概述 (Summary)

产品名称 (Production Name)	PROKR2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Mouse
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PROKR2
别名 (Alternative Names)	HH3; KAL3; PKR2; GPRg2; GPR73b; GPR73L1; dJ680N4.3
基因 ID (Gene ID)	128674.0
蛋白 ID (SwissProt ID)	Q8NFJ6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Observed MW: 42 kDa

研究背景 (Background)

prokineticin receptor 2(PROKR2) Homo sapiens Prokineticins are secreted proteins that can promote angiogenesis and induce

产品名称: PROKR2 Rabbit Polyclonal Antibody
产品货号: AP Rab20429




strong gastrointestinal smooth muscle contraction. The protein encoded by this gene is an integral membrane protein and G

protein-coupled receptor for prokineticins. The encoded protein is similar in sequence to GPR73, another G protein-coupled receptor for prokineticins. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .