

产品名称: METRL Rabbit Polyclonal Antibody
产品货号: APRab20440



产品概述 (Summary)

产品名称 (Production Name)	METRL Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide.
储存溶液 (Buffer)	sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	METRNL
别名 (Alternative Names)	-
基因 ID (Gene ID)	284207.0
蛋白 ID (SwissProt ID)	Q641Q3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	Calculated MW: 34 kDa


研究背景 (Background)

产品名称: METRL Rabbit Polyclonal Antibody
产品货号: APRab20440




研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4°over night



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2,

Tris-

EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

注意事项 (Note)

For research use only .