

产品名称: NRAM2 Rabbit Polyclonal Antibody
产品货号: APRab20444



产品概述 (Summary)

产品名称 (Production Name)	NRAM2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NRAM2
别名 (Alternative Names)	SLC11A2 DCT1 DMT1 NRAMP2 OK/SW-cl.20
基因 ID (Gene ID)	4891.0
蛋白 ID (SwissProt ID)	P49281.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Observed MW: 60 kDa

研究背景 (Background)

This gene encodes a member of the solute carrier family 11 protein family. The product of this gene transports divalent


产品名称: NRAM2 Rabbit Polyclonal Antibody
产品货号: APRab20444



metals and is involved in iron absorption. Mutations in this gene are associated with hypochromic microcytic anemia with iron overload. A related solute carrier family 11 protein gene is located on chromosome 2. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Apr 2010],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .