

产品名称: RIP (Phospho Ser166) Rabbit Polyclonal Antibody
产品货号: APRab20453

产品概述 (Summary)

产品名称 (Production Name)	RIP (Phospho Ser166) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RIPK1
别名 (Alternative Names)	RIPK1;RIP;RIP1Receptor-interacting protein 1 ;RIP-1 ;
基因 ID (Gene ID)	8737.0
蛋白 ID (SwissProt ID)	Q13546.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	Observed MW: 73 kDa

研究背景 (Background)

Essential adapter molecule for the activation of NF-kappa-B. Following different upstream signals (binding of inflammatory


产品名称: RIP (Phospho Ser166) Rabbit Polyclonal Antibody
产品货号: APRab20453



cytokines, stimulation of pathogen recognition receptors, or DNA damage), particular RIPK1-containing complexes are formed, initiating a limited number of cellular responses.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human cervical carcinoma. 1, Antibody was diluted at 1:200(4° overnight).

2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).

注意事项 (Note)

For research use only .