

产品名称: FUNDC1 Rabbit Polyclonal Antibody
产品货号: APRab20477



产品概述 (Summary)

产品名称 (Production Name)	FUNDC1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02%
储存溶液 (Buffer)	sodium azide.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FUNDC1
别名 (Alternative Names)	FUNDC1 ;FUN14 domain-containing protein 1
基因 ID (Gene ID)	139341.0
蛋白 ID (SwissProt ID)	Q8IVP5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:300, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Observed MW: 17 kDa

研究背景 (Background)

产品名称: FUNDC1 Rabbit Polyclonal Antibody
产品货号: APRab20477




This gene encodes a protein with a FUN14 superfamily domain. The function of the encoded protein is not known.

[provided


by RefSeq, Sep 2011],

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of AD293 cells using FUNDC1 Polyclonal Antibody diluted at 1:500.



Immunofluorescence analysis of A549. 1, primary Antibody(red) was diluted at 1:200(4°C overnight). 2, Goat Anti Rabbit IgG (H&L) -Alexa Fluor 594 Secondary antibody was diluted at 1:1000(room temperature, 50min).3, DAPI(blue) 10min.

注意事项 (Note)

For research use only .