

产品名称: MCM6 Rabbit Monoclonal Antibody
产品货号: AMRe01410



产品概述 (Summary)

产品名称 (Production Name)	MCM6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MCM6
别名 (Alternative Names)	minichromosome maintenance complex component 6; Mis5; P105MCM; MCG40308
基因 ID (Gene ID)	4175
蛋白 ID (SwissProt ID)	Q14566.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 93 kDa; Observed MW: 93 kDa

研究背景 (Background)

产品名称: MCM6 Rabbit Monoclonal Antibody
产品货号: AMRe01410




Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of MCM6 in Jurkat lysates using MCM6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of MCM6 (green) in Hela using MCM6 antibody, and DAPI (blue)

产品名称: MCM6 Rabbit Monoclonal Antibody
产品货号: AMRe01410



注意事项 (Note)

For research use only .