

产品名称: Rab8A Rabbit Monoclonal Antibody
产品货号: AMRe01423



产品概述 (Summary)

产品名称 (Production Name)	Rab8A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB8A
别名 (Alternative Names)	RAB8A; member RAS oncogene family; MEL; RAB8
基因 ID (Gene ID)	4218
蛋白 ID (SwissProt ID)	P61006.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

研究背景 (Background)

The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to

产品名称: Rab8A Rabbit Monoclonal Antibody
产品货号: AMRe01423




their fusion with membranes.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of RAB8A in K562, C6, 3T3, HeLa lysates using Rab8A antibody.



Immunocytochemistry analysis of Rab8A(green) in HeLa using Rab8A antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using RAB8A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Rab8A Rabbit Monoclonal Antibody
产品货号: AMRe01423



注意事项 (Note)

For research use only .