

产品名称: MSP Rabbit Monoclonal Antibody
产品货号: AMRe01440



产品概述 (Summary)

产品名称 (Production Name)	MSP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STK4
别名 (Alternative Names)	serine threonine kinase 3/4; KRS1; MST2/KRS2; MST1; YSK3; TIAC
基因 ID (Gene ID)	6789
蛋白 ID (SwissProt ID)	Q13043.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000, ICC/IF 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 60 kDa

研究背景 (Background)

Stress-activated, pro-apoptotic kinase which, following caspase-cleavage, enters the nucleus and induces chromatin

产品名称: MSP Rabbit Monoclonal Antibody
产品货号: AMRe01440




condensation followed by internucleosomal DNA fragmentation.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunocytochemistry analysis of MSP (green) in K562 using MSP antibody, and DAPI (blue).



Western blot analysis of MSP/MST1 in K562, rat Brain, C6, 3T3, HeLa lysates using MSP/MST1 antibody.

注意事项 (Note)

For research use only .