

产品名称: Annexin VI Rabbit Monoclonal Antibody
产品货号: AMRe01445



产品概述 (Summary)

产品名称 (Production Name)	Annexin VI Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA6
别名 (Alternative Names)	p68; p70; ANX6; CBP68; CPB-II
基因 ID (Gene ID)	309
蛋白 ID (SwissProt ID)	P08133.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 76 kDa; Observed MW: 70 kDa

研究背景 (Background)

May associate with CD21. May regulate the release of Ca²⁺ from intracellular stores.MiscellaneousSeems to bind one

产品名称: Annexin VI Rabbit Monoclonal Antibody
产品货号: AMRe01445




calcium ion with high affinity.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of Annexin VI in C6 lysates using Annexin VI antibody.



Immunocytochemistry analysis of Annexin VI(green) in HeLa using Annexin VI antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using Annexin VI antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Annexin VI Rabbit Monoclonal Antibody
产品货号: AMRe01445



注意事项 (Note)

For research use only .