

产品名称: TRMT2A Rabbit Monoclonal Antibody
产品货号: AMRe01446



产品概述 (Summary)

产品名称 (Production Name)	TRMT2A Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRMT2A
别名 (Alternative Names)	Hpa11 tiny fragments locus 9c protein
基因 ID (Gene ID)	27037
蛋白 ID (SwissProt ID)	Q8IZ69.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 69 kDa; Observed MW: 69 kDa

研究背景 (Background)

S-adenosyl-L-methionine-dependent methyltransferase that catalyzes the formation of 5-methyl-uridine in tRNAs and

产品名称: TRMT2A Rabbit Monoclonal Antibody
产品货号: AMRe01446




some mRNAs.


研究领域 (Research Area)

Epigenetics and Nuclear Signaling


图片 (Image Data)



Western blot analysis of HTF9C/TRMT2A in K562, HeLa lysates using TRMT2A antibody.



Immunocytochemistry analysis of HTF9C/TRMT2A(green) in HeLa using HTF9C/TRMT2A antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using HTF9C antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: TRMT2A Rabbit Monoclonal Antibody
产品货号: AMRe01446



注意事项 (Note)

For research use only .